



EMBEDDED DICOM WEB VIEWER



Perennity

ACCESSBOX




Review, Transfer, Share & Protect your Medical Images

● APPLIANCE ● EMBEDDED ● ENTERPRISE ● WINDOWS/MAC CLIENT

All-in-One for Medical Imaging


Transfer images between institutions

Securely exchange studies between Imaging & Diagnostic centers, clinics, hospitals and radiologists. No VPN or cloud services required.



Local Storage/Archive

Use the AccessBox as a genuine DICOM server to store, query & retrieve your DICOM images and reports.



Review studies using any WEB browser

The AccessBox embeds a full-featured DICOM WEB viewer. Review your images using Windows, MAC, Linux, iOS and Android.



Protect your images

Replicate your studies and extend the storage capacity of the AccessBox using any local NAS or any cloud repository.



Benefits

Institutions

- Easily share studies between hospitals, clinics and imaging & diagnostic centers
- Dispatch studies to Radiologists outside of the facility for diagnostic purpose
- Archive images to secure internal/external/cloud-based Data Center
- Works with limited bandwidth and even unstable Internet connection

Radiologists

- Quickly access studies from anywhere, anytime
- Easily collect studies from imaging center for diagnostic
- Send report using any application
- Publish studies to Web Portal for Referring Physicians and Patients
- No IT expertise needed
- Available as software portable client for Window & Mac

Patients

- Faster diagnostics
- Access to reports & images
- Download original images for secondary diagnostics
- Economical

Physicians

- Get instant access to images and reports of their patients
- Email notification when study is available for review
- Securely access images from any computer, tablet or smartphone

Four versions to meet your needs

Standalone



Out-of-the-box appliance
Preconfigured
1TB integrated hard drive
Extendable using any NAS
Up to 5 nodes

Embedded



Hosted on Synology NAS
Embedded RAID
High storage capacity
Highly expandable
Up to 20 nodes

Enterprise




Hosted on MS Windows
Server 2012/2016
iMAGE Portal included
Unlimited storage capacity
Unlimited number of nodes

Client



Software plug-in
Windows or Mac
Send/Receive images with
your workstation to/from
any AccessBox

The AccessBox optimizes Patient Data Exchange and Archiving





Perennity

ACCESSBOX

All-in-One for Medical Imaging

SECURELY REVIEW, TRANSFER, SHARE & PROTECT

HARBTEK
SYSTEMS

Distributed By:

Harb-TEK Systems
info@harb-tek.com www.harb-tek.com



Perennity

EMEA

Brussels (Belgium)
Tel: +32 (0)2/421.59.60
Fax: +32 (0)2/421.59.69
<http://www.perennitysoft.eu>
sales@perennitysoft.eu

Perennity

AMERICAS

Miramar (Florida - USA)
Tel: +1 954-734-2836
Fax: +1 954-874-1699
<http://www.perennitysoft.com>
sales@perennitysoft.com

Perennity

ASIA - PACIFIC

Singapore
Phone: +65 67432325
Fax: +65 67437032
<http://www.perennitysoft.com>
asiapac@perennitysoft.com